

# **USER GUIDE**

# Topbar & Clinical Desktop Version 1.0

Issued: 02.08.2017

Author: Damian Cutler & Jon Rolph



## Contents

Clinical Desktop / Starting Topbar	3
Topbar	7
Topbar applications	7
Adding and removing apps	7
Data Cleansing app	9
Demographic tab	9
Clinical tab	.10
Indications tab	.11
Filters tab	.12
Patient Health Summary app	.12
MBS Items app	.13
Guidelines tab	.13
Activities tab	.14
Claims tab	.14
CAT Prompts app	.15
Notifications tab	.15
Settings tab	.16
Waiting Room app	.16
Creating CAT Plus prompts	.18



## Clinical Desktop / Starting Topbar


1. Logon to the clinical desktop via the shortcut



2. The clinical desktop is where all your clincal applications are located. It is made to look different to your local desktop so you are never in doubt about where you are.




3. Topbar should automatically download when you login.




4. Topbar will then install/update itself under your profile

2/08/2017 Page: 3 Version 1.0





5. Accept the Topbar disclaimer.



2/08/2017 Page: 4 Version 1.0



6. You will then see the Topbar release notes




7. Topbar will then open at the top of the screen, please use your Pen CS (used for CAT) login credentials to login. Contact your supervisor if you require credentials.



2/08/2017 Page: 5 Version 1.0



8. Once you are logged into Topbar, open Medical Director Clinical. When you open a patient's record Topbar will indicate any actions or prompts that it finds for that patient.



2/08/2017 Page: 6 Version 1.0



### **Topbar**

#### **Topbar applications**

Topbar currently has the following apps:

**Data Cleansing App Overview** 

Patient Health Summary PHS App Overview

**MBS App Overview** 

**CAT Prompts App Overview** 


**Waiting Room App Overview** 

These apps are displayed in Topbar only if the currently open patient has any actions recommended. Each app can be opened manually by changing the settings in the dashboard as described here (or below): Topbar Dashboard

#### Adding and removing apps

You can configure your Topbar to display the apps you want, in the order you want.


1. Click on the Topbar icon in the far left of the Topbar, and select edit app.



2. You can then disable apps by hitting on the cross on the top right of each app, they will then appear in the Disabled Applications pane on the righthand side, you can pin an app and this app will then always display in the Topbar or you can change the order of the apps by left click, hold and drag the apps into the order you would like.

2/08/2017 Page: 7 Version 1.0

# Topbar & Clinical Desktop ITERAMS



2/08/2017 Page: 8 Version 1.0



#### Data Cleansing app

The data cleansing app has four tabs:

data and also

• Demographic - Displays any missing demographical data for the open patient

• Clinical - Displays any missing clinical / accreditation data for the open patient

• Indications - Indicates if a patient potentially has a chronic disease

• Filters - Where you can turn off demographic/clinical/indicated fields for

**Topbar** 

#### Demographic tab

This tab willI show you any missing patient demographics on the open patient.


You can either:

Add in Clinical System - Will open the location on the patient record to

enter the data

Defer
 Will offer you a future date to defer entering the

an optional reason for defering



2/08/2017 Page: 9 Version 1.0



#### Clinical tab

This tab will show you any missing clinical / accreditation items on the open patient.

#### You can either:


• Add in Clinical System - Will open the location on the patient record to

enter the data

Defer

 Will offer you a future date to defer entering the data and also

an optional reason for defering



2/08/2017 Page: 10 Version 1.0

#### Indications tab

This tab will show you potential chronic diseases the patient may have for which there is no active diagnosis.

#### You can either:


- Add in Clinical System
- Defer entering the data
- Confirm indicated condition not exist and does not exist
- Topbar or

- Will open the location on the patient record to enter the data
- Will offer you a future date to defer to before

and also an optional reason for defering
Clinician can confirm that indicated condition does

patient will no longer be flagged for the condition in

CAT4 reporting. Reviwed conditioner will appear under review at the bottom of the tab and can be reinstated.




2/08/2017 Page: 11 Version 1.0



#### Filters tab


This tab will let you toggle which items you want active on the Demographics, Clinical and Indicated tabs.

You can turn each item off and also have the option with the Clinic items to set a date when these items were last updated.



#### Patient Health Summary app

Gives an overview of the patient record starting with most current items on the left and historic items to the right which can be scrolled across using the bottom scroll bar. No input on this tab.



2/08/2017 Page: 12 Version 1.0



#### MBS Items app

The MBS has 3 tabs:

 Guidelines - Shows relevant MBS items for the patient and if they are up to date

 Reccommended - Reccommends activities for the patient based upon Medicare guidelines

• Claims - Shows MBS items Ready for claiming, items Claimable if

activities are completed, Up-to-date MBS items and items


Not eligible for claiming

#### Guidelines tab

This tab will show relevant and not relevant care items for the patient.

Once an item is selected, you can either:

- Visit MBS Online for further information about the item
- Add in the clinical system




2/08/2017 Page: 13 Version 1.0




#### Activities tab

This tab will show recommended activities and the claim eligibility (<u>based on our records only</u>) for a patient and allows a clinician to add the activity into the clinical system.



#### Claims tab

Shows MBS items Ready for claiming, items Claimable if activities are completed, Up-to-date MBS items and items Not eligible for claiming. Claimable gives you the option to add the activity to the clinical system to claim.



2/08/2017 Page: 14 Version 1.0



#### **CAT Prompts app**

This app allows a clinician to create their own triggers for measures against patient records. These can be added through the Daily View of the CAT4 and any patient meeting the criteria, will be flagged.


#### Notifications tab

This will flag against any patient that meets the criteria of a prompt.

#### You can either:

- Dismiss Will dismiss the prompt against the open patient
- Defer Will offer you a future date to defer entering the data and also an optional reason for defering

Any any Deffered or Dismissed prompts will display under the respective headings.




2/08/2017 Page: 15 Version 1.0



#### Settings tab


The settings tab allows you to view details of selected prompts and also toggle whether you want them displayed or not.



#### Waiting Room app


This app flags any patients missing Demographic or Clinical items in the same way as the tabs from the Data Cleansing app except it will flag any patient put into the PracSoft waiting room instead of the open patient.

Same options apply as the Data Cleansing app and the Waiting Room app also has a filters tab for toggling filters



2/08/2017 Page: 16 Version 1.0

# Topbar & Clinical Desktop IT RAMS



2/08/2017 Page: 17 Version 1.0



## **Creating CAT Plus prompts**

You can create your own patient prompts from the CAT4 but only from the Daily View presently.

1. Login to the CAT4. If you don't have an account, please talk with your Supervisor.



#### 2. Select Daily View




For use of the CAT4, please see <a href="here">here</a> for an IT@RAMS getting started guide or additional resources from Pen CS can be found <a href="here">here</a> for CAT Plus & <a href="here">here</a> for the CAT4.


2/08/2017 Page: 18 Version 1.0



3. Once you have added your filters to create the criteria you wish patients to be prompted against, hit the export button.




4. Use the list box at the bottom of the patient Reidentification window to select "Prompt at Consult – Topbar" and hit the Go button.



2/08/2017 Page: 19 Version 1.0



5. You will then need to add a meaningful name to your prompt.



2/08/2017 Page: 20 Version 1.0